

CARE FLOW

Transurethral *Soft* 100% silicone catheter

2-way balloon catheter with nelaton tip

Product specifications:

- ✓ In CH sizes 06 – 10 the catheter comes with a stabilising guiderail.
- ✓ Very high wearing comfort, due to the soft-silicone material of the catheter shaft.
- ✓ Safe fixation of the catheter in the bladder due to the balloon.
- ✓ Use our blocking syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The X-ray line allows for radiological imaging of the catheter.
- ✓ The catheter can be used for long-term drainage.



Nelaton tip

Size	Balloon	Length	Article number	PZN
CH 06	03 ml	31 cm	01/01/01/06/03/GHC	17883185
CH 08	05 ml	31 cm	01/01/01/08/05/GHC	17883216
CH 10	05 ml	31 cm	01/01/01/10/05/GHC	17883222
CH 12	10 ml	40 cm	01/01/01/12/10/GHC	17883239
CH 14	10 ml	40 cm	01/01/01/14/10/GHC	17883251
CH 16	30 ml	40 cm	01/01/01/16/30/GHC	17883268
CH 18	30 ml	40 cm	01/01/01/18/30/GHC	17883274
CH 20	30 ml	40 cm	01/01/01/20/30/GHC	17883280
CH 22	30 ml	40 cm	01/01/01/22/30/GHC	17883297
CH 24	30 ml	40 cm	01/01/01/24/30/GHC	17883305
CH 26	30 ml	40 cm	01/01/01/26/30/GHC	17883311

Packaging Unit: 10 pieces
Packaging: sterile



Balloon volume: 03-30 ml

Pictures may differ from the original