



### 2-way balloon catheter

- ✓ In CH sizes 06 – 10 the catheter comes with a stabilising guiderail.
- ✓ Due to the high balloon capacity of up to 30 ml, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The X-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for long-term drainage.



Nelaton tip



Balloon volume: 03-30 ml

Size	Balloon	Length	Article number	PZN
CH 06	03 ml	27 cm	01/01/11/06/03/GHC	01150196
CH 08	05 ml	27 cm	01/01/11/08/05/GHC	01150204
CH 10	05 ml	27 cm	01/01/11/10/05/GHC	01150210
CH 12	10 ml	40 cm	01/01/11/12/10/GHC	01150368
CH 14	10 ml	40 cm	01/01/11/14/10/GHC	01150374
CH 16	30 ml	40 cm	01/01/11/16/30/GHC	01150380
CH 18	30 ml	40 cm	01/01/11/18/30/GHC	01150397
CH 20	30 ml	40 cm	01/01/11/20/30/GHC	01150428
CH 22	30 ml	40 cm	01/01/11/22/30/GHC	01150434
CH 24	30 ml	40 cm	01/01/11/24/30/GHC	01150440
CH 26	30 ml	40 cm	01/01/11/26/30/GHC	01150457

**Packaging Unit:** 10 pieces

**Packaging:** sterile