


CARE FLOW

Transurethral 100% silicone catheter

2-way balloon catheter with Tiemann tip

Product specifications:

- ✓ An atraumatically dilating ball tip ensures a smooth passage of the prostatic fossa.
- ✓ The catheter offers high patient comfort.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The X-ray line allows for radiological imaging of the catheter.
- ✓ The catheter can be used for long-term drainage.



Tiemann tip

| Size | balloon | length | Article number | PZN |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml | 40 cm | 01/01/12/12/05/GHC | 01150606 |
| CH 14 | 20 ml | 40 cm | 01/01/12/14/20/GHC | 01150629 |
| CH 16 | 30 ml | 40 cm | 01/01/12/16/30/GHC | 01150635 |
| CH 18 | 30 ml | 40 cm | 01/01/12/18/30/GHC | 01150641 |
| CH 20 | 30 ml | 40 cm | 01/01/12/20/30/GHC | 01150670 |
| CH 22 | 30 ml | 40 cm | 01/01/12/22/30/GHC | 01150687 |
| CH 24 | 30 ml | 40 cm | 01/01/12/24/30/GHC | 01150693 |

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 5-30 ml

Pictures may differ from the original