

CARE FLOW

Transurethral siliconised latex catheter

2-way balloon catheter with nelaton tip

Product specifications:

- ✓ Made of high-quality, siliconized medical latex.
- ✓ The catheter offers high patient comfort.
- ✓ Due to the high balloon capacity, the catheter is safety fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The catheter is suited for short-term drainage of up to 5 days only.
- ✓ This product contains natural latex which may cause allergic reactions.



Nelaton tip

Size	Balloon	Length	Article number	PZN
CH 12	10 ml	40 cm	01/01/13/12/10/GHC	01150701
CH 14	10 ml	40 cm	01/01/13/14/10/GHC	01150859
CH 16	30 ml	40 cm	01/01/13/16/30/GHC	01150865
CH 18	30 ml	40 cm	01/01/13/18/30/GHC	01150871
CH 20	30 ml	40 cm	01/01/13/20/30/GHC	01150888
CH 22	30 ml	40 cm	01/01/13/22/30/GHC	01150919

Packaging Unit: 10 pieces
Packaging: sterile



Balloon volume: 10-30 ml

Pictures may differ from the original