


CARE FLOW

Transurethral 100% silicone irrigation catheter

3-way balloon catheter with mercier tip

Product specifications:

- ✓ An atraumatic and kink resistant Mercier tip ensures a smooth passage of the prostatic fossa.
- ✓ Optimal drainage and irrigation due to two staggered drainage eyes and an open tip.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The x-ray tip allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for post-operative long-term drainage.



Mercier tip

Size	Balloon	Length	Article number	PZN
CH 18	40 ml	40 cm	01/01/18/18/40/GHC	01151511
CH 20	50 ml	40 cm	01/01/18/20/50/GHC	01151528
CH 22	60 ml	40 cm	01/01/18/22/60/GHC	01151534
CH 24	70 ml	40 cm	01/01/18/24/70/GHC	01151540

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 40-70 ml

Pictures may differ from the original