


CARE FLOW

Transurethral PVC single-use male catheter

Transparent with funnel, nelaton and tiemann tip

Product specifications:

- ✓ Disposable female catheter with two staggered atraumatic drainage eyes.
- ✓ Rounded and closed Nelaton or Tiemann tip, for a smooth passage of the prostatic fossa.
- ✓ Soft PVC material.
- ✓ The catheter is for single use only.



Nelaton tip



Tiemann tip

Nelaton tip

Size	Length	Article number	PZN
CH 08	40 cm	01/01/20/08/00/GHC	14289493
CH 10	40 cm	01/01/20/10/00/GHC	08871349
CH 12	40 cm	01/01/20/12/00/GHC	08871355
CH 14	40 cm	01/01/20/14/00/GHC	08871361
CH 16	40 cm	01/01/20/16/00/GHC	08871378
CH 18	40 cm	01/01/20/18/00/GHC	08871384
CH 20	40 cm	01/01/20/20/00/GHC	08871390

Tiemann tip

Size	Length	Article number	PZN
CH 10	40 cm	01/01/21/10/00/GHC	01077312
CH 12	40 cm	01/01/21/12/00/GHC	01929620
CH 14	40 cm	01/01/21/14/00/GHC	01929637
CH 16	40 cm	01/01/21/16/00/GHC	01929666
CH 18	40 cm	01/01/21/18/00/GHC	01971792
CH 20	40 cm	01/01/21/20/00/GHC	09219740

Packaging Unit: CH 08-18: 50 pieces each, CH 20: 30 pieces

Packaging: sterile

Pictures may differ from the original