

CARE FLOW

Transurethral 100% silicone irrigation catheter

3-way balloon catheter with nelaton tip

Product specifications:

- ✓ An atraumatic Nelaton tip ensures a smooth passage of the prostatic fossa.
- ✓ Optimal drainage and irrigation performance, due to two staggered drainage eyes.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for post-operative long-term drainage.



Nelaton tip

Size	Balloon	Length	Article number	PZN
CH 18	30 ml	40 cm	01/01/22/18/30/GHC	15618237
CH 20	30 ml	40 cm	01/01/22/20/30/GHC	15618243
CH 22	30 ml	40 cm	01/01/22/22/30/GHC	15618266
CH 24	30 ml	40 cm	01/01/22/24/30/GHC	15618295

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 30 ml

Pictures may differ from the original