

CARE FLOW

Transurethral PVC single-use female catheter

Transparent with funnel and nelaton tip, DEHP-free

Product specifications:

- ✓ Disposable female catheter with two staggered atraumatic drainage eyes.
- ✓ Rounded and closed tip.
- ✓ Soft PVC material.
- ✓ DEHP-free.
- ✓ The catheter is for single use only.



Nelaton tip

Size	Length	Article number	PZN
CH 06	20 cm	01/01/35/06/00/GHC	18320969
CH 08	20 cm	01/01/35/08/00/GHC	18320975
CH 10	20 cm	01/01/35/10/00/GHC	18320981
CH 12	20 cm	01/01/35/12/00/GHC	18320998
CH 14	20 cm	01/01/35/14/00/GHC	18321006
CH 16	20 cm	01/01/35/16/00/GHC	18321035

Packaging Unit: 100 pieces

Packaging: sterile



Laterally displaced eyes

Abbildungen können vom Original abweichen à Pictures may differ from the original