

Additional equipment

Accessories for the catheter replacement

Next to our assortment covering catheters as well as drainage- and holding- systems, we also provide necessary accessories, to maintain a satisfactory supply of revulsive incontinence aids and ensuring an optimal replacement of a catheter, which all helps to maintain a high level of supply quality.

We want to highlight:

- The **blocking syringe containing a 10% glycerine solution, 10 ml** (Art.-Nr. 01/01/99/KB/10/GHC), ensures a steady volume in the catheter balloon during the wearing period, by preventing the liquid medium from diffusing out of the balloon. Thereby the balloon catheter stays fixated, and the risk of a dislocated catheter can almost be excluded.
- The **lubricant LUBRICARE, 6 ml** (Art.-Nr. 01/01/99/GG/001), adhere excellently onto the mucous membrane and has excellent lubrication properties, which ensures a smooth replacement of a catheter.

The below illustrated products are only an extract of our comprehensive assortment. We gladly provide more information upon request.



Blocker syringe with 10% glycerine, 10 ml



Lubricant LUBRICARE, 06 ml



Catheter plug



Octenisept, 15 ml

Article	Article number	PZN	VE
1. Blocker syringe	01/01/99/KB/10/GHC	06718980	1VE = 200 pieces
2. LUBRICARE Lubricant	01/01/99/GG/06/001	13360220	1VE = 100 pieces
LUBRICARE Lubricant	01/01/99/GG/06/25	16938435	1VE = 4 x 25 pieces
Endosgel Lubricant	01/01/99/HG/EG/GHC	02437623	1VE = 250 pieces
3. Catheter plug	01/01/99/ZB/ST/GHC	14354919	1VE = 20 pieces
4. Octenisept	01/01/99/OCT/15/GHC	04804008	1VE = 100 pieces

Packaging: sterile

Pictures may differ from the original