

CARE SET

Catheter insertion kit Comfort Eco

Product specifications:

- ✓ The Comfort Eco Set comprises an inner and an outer package.
- ✓ The outer package of the basic set contains: one sterile 10 ml syringe for unblocking the inserted catheter, mucous membrane disinfectant Octenisept 15 ml, lubricant LUBRICARE 6 ml.
- ✓ The items in the inner, sterilised package are used to place a new transurethral catheter.



Article	Article number	PZN
Comfort Eco	01/01/99/KL/CE/GHC	16134211

Components of the catheter insertion kit are:

1 Tray 	6 Swabs 	1 Pair of tweezers 	10 ml Syringe 	Mucous membrane 
1 Blocker syringe with 10% glycerine 	3 Nitrile gloves 	1 Fenestrated drape 	2 Underpads 	Octenisept, 15 ml Lubricant 

Packaging Unit: 35 pieces
Packaging: sterile

Pictures may differ from the original