

# CARE SET

## Catheter insertion kit Comfort and Comfort Plus

### Product specifications:

- ✓ The Catheter insertion Kits comprises an inner and an outer package.
- ✓ The outer package of the basic set contains: one sterile 10 ml syringe for unblocking the inserted catheter a mucous membrane disinfectant Octenisept 15 ml.
- ✓ Additional included in Comfort Plus: lubricant Endosgel 6 ml.
- ✓ The items in the inner, sterilised package are used to place a new transurethral catheter.



Article	Article number	PZN
Comfort	01/01/99/KL/CO/GHC	14209815
Comfort Plus	01/01/99/KL/CP/GHC	14209821

### Components of the catheter insertion kit are:

1 Tray



1 Blocker syringe with 10% glycerine




6 Swabs



3 Nitrile gloves



1 Pair of tweezers



1 Fenestrated drape



10 ml Syringe



2 Underpads



Mucous membrane disinfectant



Octenisept 15 ml

Additional included in Comfort Plus:



Lubricant Endosgel 6 ml

**Packaging Unit:** 35 pieces

**Packaging:** sterile

Pictures may differ from the original