

CARE SET

Catheter insertion kit Comfort and Comfort Plus

Product specifications:

- ✓ The Catheter insertion Sets comprises an inner and an outer package.
- ✓ The outer package of the basic set contains: one sterile 10 ml syringe for unblocking the inserted catheter, mucous membrane disinfectant Octenisept 15 ml.
- ✓ additional included in Comfort Plus: lubricant Endosgel 6 ml.
- ✓ The items in the inner, sterilised package are used to place a new transurethral catheter.



Article	Article number	PZN
Comfort	01/01/99/KL/CO/GHC	14209815
Comfort Plus	01/01/99/KL/CP/GHC	14209821

Components of the catheter insertion kit are:


1 tray



6 swabs



1 pair of tweezers



10 ml syringe



mucous membrane disinfectant



Octenisept 15 ml

additional included in Comfort Plus:



lubricant Endosgel 6 ml

1 blocker syringe with 10% glycerine



3 nitrile gloves



1 fenestrated drape



2 underpads



Packaging Unit: 35 pieces

Packaging: sterile

Pictures may differ from the original