


CARE FLOW

Suprapubic 100% silicone catheter *integral*

2-way balloon catheter with plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers an excellent wearing comfort.
- ✓ We recommend using our blocking syringe containing a 10% glycerine solution, for a safe fixation of the catheter
- ✓ 100% silicone for excellent biocompatibility and high resistance to encrustation.
- ✓ Short and central open tip.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.



Central open tip, integral balloon

25

Size	Balloon	Length	Article number	PZN
CH 12	05 ml	41 cm	01/02/09/12/05/GHC	16840937
CH 14	05 ml	41 cm	01/02/09/14/05/GHC	16840943
CH 16	10 ml	41 cm	01/02/09/16/10/GHC	16840972
CH 18	10 ml	41 cm	01/02/09/18/10/GHC	16840989
CH 20	10 ml	41 cm	01/02/09/20/10/GHC	16840995
CH 22	10 ml	41 cm	01/02/09/22/10/GHC	16841003

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original