

# CARE FLOW

## Suprapubic 100% silicone catheter *integral blue*

### 2-way balloon catheter with plug

#### Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers an excellent wearing comfort.
- ✓ We recommend using our blocking syringe containing a 10% glycerine solution, for a safe fixation of the catheter.
- ✓ 100% silicone for excellent biocompatibility and high resistance to encrustation.
- ✓ Short and central open tip.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optimal placement.



*Central open, integral balloon*

Size	Balloon	Length	Article number	PZN
CH 12	05 ml	41 cm	01/02/11/12/05/GHC	16841026
CH 14	05 ml	41 cm	01/02/11/14/05/GHC	16841032
CH 16	10 ml	41 cm	01/02/11/16/10/GHC	16841049
CH 18	10 ml	41 cm	01/02/11/18/10/GHC	16841061
CH 20	10 ml	41 cm	01/02/11/20/10/GHC	16841078
CH 22	10 ml	41 cm	01/02/11/22/10/GHC	16841084

**Packaging Unit:** 10 pieces

**Packaging:** sterile



*Balloon volume: 05-10 ml*

Pictures may differ from the original