

CARE FLOW UNIBAL

Suprapubic *integral* 100% silicone catheter

2-way integral balloon catheter

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers high patient comfort.
- ✓ High drainage performance, due to the high lumen and the central open tip.
- ✓ We recommend using or blocking syringe containing a 10% glycerine solution, for a safe fixation of the catheter
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Short and central open tip.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.



Central open tip, integral balloon

Size	Balloon	Length	Article number	PZN
CH 10	03 ml	40 cm	01/02/12/10/03/GHC	01152568
CH 12	05 ml	40 cm	01/02/12/12/05/GHC	01152574
CH 14	05 ml	40 cm	01/02/12/14/05/GHC	01152605
CH 16	10 ml	40 cm	01/02/12/16/10/GHC	01152611
CH 18	10 ml	40 cm	01/02/12/18/10/GHC	01152628
CH 20	10 ml	40 cm	01/02/12/20/10/GHC	01152700
CH 22	10 ml	40 cm	01/02/12/22/10/GHC	01152717

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 03-10 ml

Pictures may differ from the original