

CARE FLOW UNIBAL

Suprapubic change set, incl. guidewire


2-way integral balloon catheter

Product specifications:

This charge set contains:

- 100% silicone *integral* balloon catheter
- Teflon-coated guidewire

- ✓ Integral balloon catheter with seamlessly incorporated balloon into the shaft.
- ✓ High drainage performance, due to the high lumen and the central open tip.
- ✓ Stable teflon-coated guidewire.
- ✓ Safe fixation of the catheter, due to the integral balloon.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ Teflon-coated guidewire, length 80 cm, 3 cm straight and flexible tip, thickness 0,035".



Central open tip, integral balloon

Guidewire

| Wire | Tip | Length | Article number | PZN |
|--------|------------|--------|--------------------|----------|
| coated | 3 cm flex. | 80 cm | 01/02/18/80/35/GHC | 14166431 |

Packaging Unit: 1 piece

Packaging: sterile

Change set, incl. Guidewire

| Size | Balloon | Length | Article number | PZN |
|-------|---------|--------|--------------------|----------|
| CH 10 | 03 ml | 40 cm | 01/02/16/10/03/GHC | 01153154 |
| CH 12 | 05 ml | 40 cm | 01/02/16/12/05/GHC | 01153160 |
| CH 14 | 05 ml | 40 cm | 01/02/16/14/05/GHC | 01153183 |

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 03-05 ml

Pictures may differ from the original