


BLUECATH

Suprapubic *integral* 100% silicone catheter

2-way integral balloon catheter with dilatation tip and plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ Safe fixation of the catheter, due to the integral balloon.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Dilatation tip with central bore.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optional placement.



Dilatation tip with central bore, integral balloon

27

Size	Balloon	Length	Article number	PZN
CH 12	05 ml	40 cm	01/02/21/12/05/GHC	11850948
CH 14	05 ml	40 cm	01/02/21/14/05/GHC	11850954
CH 16	10 ml	40 cm	01/02/21/16/10/GHC	11850983
CH 18	10 ml	40 cm	01/02/21/18/10/GHC	11851008
CH 20	10 ml	40 cm	01/02/21/20/10/GHC	11851014
CH 22	10 ml	40 cm	01/02/21/22/10/GHC	11851020

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original