

# BLUECATH

## Suprapubic 100% silicone catheter

2-way integral balloon catheter with dilatation tip and plug

### Product specifications:

- ✓ The catheter is equipped with a specially formed conical tip preventing urethral stricture occurring with conventional catheters. Consequently, there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers high patient comfort.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.



*Dilatation tip with central bore, integral balloon*

Size	balloon	length	Article number	PZN
CH 12	05 ml	40 cm	01/02/21/12/05/GHC	11850948
CH 14	05 ml	40 cm	01/02/21/14/05/GHC	11850954
CH 16	10 ml	40 cm	01/02/21/16/10/GHC	11850983
CH 18	10 ml	40 cm	01/02/21/18/10/GHC	11851008
CH 20	10 ml	40 cm	01/02/21/20/10/GHC	11851014
CH 22	10 ml	40 cm	01/02/21/22/10/GHC	11851020

**Packaging Unit:** 10 pieces  
**Packaging:** sterile



*Balloon volume: 05-10 ml*

Pictures may differ from the original