


BLUECATH

Suprapubic 100% silicone catheter

2-way integral balloon catheter with plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon preventing urethral stricture occurring with conventional catheters. Consequently, there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers high patient comfort.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.



Central open, integral balloon

Size	balloon	length	Article number	PZN
CH 12	05 ml	40 cm	01/02/22/12/05/GHC	11851037
CH 14	05 ml	40 cm	01/02/22/14/05/GHC	11851043
CH 16	10 ml	40 cm	01/02/22/16/10/GHC	11851072
CH 18	10 ml	40 cm	01/02/22/18/10/GHC	11851089
CH 20	10 ml	40 cm	01/02/22/20/10/GHC	11851095
CH 22	10 ml	40 cm	01/02/22/22/10/GHC	11851126
CH 24	10 ml	40 cm	01/02/22/24/10/GHC	11851132

Packaging Unit: 10 pieces
Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original