


BLUECATH

Suprapubic *integral* 100% silicone catheter

2-way integral balloon catheter with plug, short

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even tickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ High drainage performance, due to the high lumen and the central open tip.
- ✓ Safe fixation of the catheter, due to the *integral* balloon.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Short, central open tip.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optional placement.



Central open tip, integral balloon

Size	Balloon	Length	Article number	PZN
CH 12	05 ml	33 cm	01/02/23/12/05/GHC	11851149
CH 14	05 ml	33 cm	01/02/23/14/05/GHC	11851155
CH 16	10 ml	33 cm	01/02/23/16/10/GHC	11851161
CH 18	10 ml	33 cm	01/02/23/18/10/GHC	11851178
CH 20	10 ml	33 cm	01/02/23/20/10/GHC	11851184
CH 22	10 ml	33 cm	01/02/23/22/10/GHC	11851190

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original