


NEPHROSTOMY

Nephrostomy silicone catheter *integral* Dilatation tip



2-way *integral* balloon catheter

- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The catheters are suitable for long-term drainage.
- ✓ Dilating nelaton tip with central bore.
- ✓ High drainage performance, due to four drainage eyes located in front of and behind the balloon.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optimal placement.
- ✓ 100% silicone.




Dilating nelaton tip with central bore, *integral* balloon

Size	Balloon	Length	Article number	PZN
CH 10	1,5 ml	40 cm	01/03/10/10/01/GHC	13851171
CH 12	2,0 ml	40 cm	01/03/10/12/02/GHC	13851188
CH 14	3,0 ml	40 cm	01/03/10/14/03/GHC	13851202
CH 16	5,0 ml	40 cm	01/03/10/16/05/GHC	13851219
CH 18	5,0 ml	40 cm	01/03/10/18/05/GHC	13851225

Packaging Unit: 10 pieces

Packaging: sterile



Balloon volume: 1,5 - 5,0 ml

Pictures may differ from the original