

CARE FLOW und CARE VENT catheter valve

For an independent and controlled emptying of the urinary bladder

Product specifications:

- ✓ The optimal solution for mobile, independent users wearing a transurethral or suprapubic catheter.
- ✓ The catheter valve continuously trains the bladder muscles so that the bladder function remains completely intact. The emptying of the urinary bladder follows the user's personal rhythm. This sustainably enhances the quality of life.
- ✓ During the night the valve can be connected to a urine bag system by means of a firmly attached silicone tube.
- ✓ By flipping the lever of the flip valve by 180° the CARE FLOW catheter valve is completely opened as well as fixed in a safe position. This prevents an unintentional closing of the valve during the night.



Arrows indicate open / closed position



Closed catheter valve



Open catheter valve, flipped by 180°

CARE FLOW catheter valve

Weight	Article number	PZN
3g	03/02/10/00/GHC	12378382

Packaging Unit: 5 pieces

Packaging: sterile

- ✓ The special construction prevents the after-dripping of urine, due to the which ensures a hygienical operation and utmost discretion.
- ✓ A colour coding of the pivot and slide mechanism makes it very easy to use single handed.

CARE VENT catheter valve

Weight	Article number	PZN
6g	03/02/10/01/GHC	01914676

Packaging Unit: 10 pieces

Packaging: non-sterile



Pictures may differ from the original