

CARE FLOW UB 750F + UB 750UF sterile / non-sterile

1-way chamber leg bag systems with flip valve for mobile users

Capacity: 750 ml

Product specifications:

- ✓ Especially users with limited hand dexterity or reduced power can safely and hygienically empty this leg bag, due to our newly developed flip valve.
- ✓ The leg bag with 50 cm tube comes with an additional stepped connector in the package. So the tube can be cut to the individual length and the leg bag can conveniently be worn on the thigh or lower leg.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.



Capacity	Tube length	Article number	PZN	Sterility
750 ml	50 cm	04/01/750F/50/GHC	12371500	sterile

Packaging Unit: 50 cm: 50 pieces

Capacity	Tube length	Article number	PZN	Sterility
750 ml	50 cm	04/01/750UF/50/GHC	12371517	non-sterile

Packaging Unit: 50 cm: 100 pieces

Packaging: each 10 pieces



Stepped connector

The special design of the connector steps provides for a reliable connection with the ring profile in the funnel end of the catheter. Needle-free urine sampling avoids needle injuries.



Doghouse

A doghouse on the bag ensures safe and hygienic positioning of the drain tap.



Fleece coating

The backside of the leg bag is coated with sweat-absorbing fleece. This effectively prevents irritation of the skin and considerably improves user comfort.



Drainage tap

The ergonomically formed flat flip valve of the drainage tube can be used easily and intuitively. To drain the leg bag the lever of the valve is flipped downwards by 180°.

Packaging: sterile / non-sterile

Fixation for leg bags on page 51-56.

Pictures may differ from the original