

# CARE FLOW UB 750 + 750US sterile / non-sterile

## 1-way chamber leg bag systems for mobile users

Capacity: 750 ml

### Product specifications:

- ✓ The leg bag with 10 cm and 30 cm tube can conveniently be worn on the thigh.
- ✓ The leg bag with 50 cm tube comes with an additional stepped connector in the package. So the tube can be cut to the individual length and the leg bag can conveniently be worn on the thigh or lower leg.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.



Capacity	Tube length	Article number	PZN	Sterility
750 ml	10 cm	04/01/750/10/GHC	06718862	sterile
750 ml	30 cm	04/01/750/30/GHC	09772046	sterile
750 ml	50 cm	04/01/750/50/GHC	06718879	sterile

**Packaging Unit:** 10 cm: 75 pieces 30 cm: 60 pieces, 50 cm: 50 pieces

Capacity	Tube length	Article number	PZN	Sterility
750 ml	10 cm	04/01/750US/10/GHC	01972716	non-sterile
750 ml	30 cm	04/01/750US/30/GHC	12341427	non-sterile
750 ml	50 cm	04/01/750US/50/GHC	00226543	non-sterile

**Packaging Unit:** 10 cm: 120 pieces 30 cm: 100 pieces, 50 cm: 100 pieces, Sub-Packaging: each 10 pieces

48



#### Stepped connector

The special design of the connector steps provides for a reliable connection with the ring profile in the funnel end of the catheter. Needle-free urine sampling avoids needle injuries.



#### Doghouse

A doghouse on the bag ensures safe and hygienic positioning of the drain tap.



#### Fleece coating

The backside of the leg bag is coated with sweat-absorbing fleece. This effectively prevents irritation of the skin and considerably improves user comfort.



#### Drain tap

The special design of the piston in the drain tap ensures safe and easy operation and prevents leakages caused by improper use of the piston.

**Packaging:** sterile / non-sterile

**Fixation for leg bags on page 51-56.**

Pictures may differ from the original