

CARE FLOW UB 2000S

Closed urine drainage systems

Filling volume: 2000 ml, soft tube

Product specifications:

- ✓ Urine drainage system with drip chamber prevents the risk of ascending germs and thus infections.
- ✓ Tube with large inner lumen made of soft material ensures very good drainage of the urine and high comfort of use.
- ✓ Twin holder with additional fastening cord enables secure attachment to different bed rails.
- ✓ Bacteria and liquid-tight ventilation filters ensure permanent pressure equalization.
- ✓ Clear scaling.
- ✓ Safe and at the same time hygienic positioning of the drain tap in the doghouse.
- ✓ Attached T-drain tap.



Volume	Tube length	Article number	PZN
2000 ml	120 cm	04/02/2000S/120/GHC	18728133
2000 ml	90 cm	04/02/2000S/90/GHC	18728179

Packaging unit: 50 pieces

Packaging: sterile

HMN: 15.25.07.0003

The **CARE FLOW BAG** serves as a sight protection for the **CARE FLOW UB 2000S** as discreet storage for the bed bags.

Product specifications:

- ✓ Handy Velcro fasteners make it quick and easy to attach the sight protection to the urine bag. A further Velcro fastener on the front allows a clean and uncomplicated usage of the drain tap.
- ✓ The sight protection is made of 100% polypropylene (PP) and is therefore resistant to organic acids, among other things.
- ✓ The sight protection is light gray, washable at 30° and therefore suitable for reuse.



Article number	PZN
04/99/01/945/GHC	09764331

Packaging unit: 10 pieces

Packaging: non-sterile

Illustrations may differ from the original