

CARE FLOW UB 2000

Closed urine drainage system

Capacity: 2000 ml

Product specifications:

- ✓ Can be fixed safely and correctly to many different types of bed frames by means of the twin holder with fixing cord.
- ✓ Bacteria- and fluid-tight venting filters ensure permanent pressure compensation and provide for optimal drainage in the system.
- ✓ Clear graduation markings allow for reliable volume control.
- ✓ A doghouse on the bag ensures safe and hygienic positioning of the drain tap.
- ✓ T-drain tap.
- ✓ Kink resistant tube.



Capacity	Tube length	Article number	PZN
2000 ml	120 cm	04/02/2000/120/GHC	06718922
2000 ml	90 cm	04/02/2000/90/GHC	09151597

Packaging Unit: 20 pieces

Packaging: sterile

The **CARE FLOW BAG** is a cover for the **CARE FLOW UB 2000** serving as a discreet way to store the urinary bedside bag.

Product specifications:

- ✓ Thanks to the convenient Velcro fastening the cover can be affixed to the urinary bag quickly and easily. A second Velcro fastening on the front allows for a clean and uncomplicated extraction and reinsertion of the drain tap.
- ✓ The cover is made from 100% polypropylene (PP) and is thus resistant against organic acids.
- ✓ The cover is light grey, washable at 30° and hence reusable.



Article number	PZN
04/99/01/945/GHC	09764331

Packaging Unit: 10 pieces

Packaging: non-sterile

Pictures may differ from the original