

CARE FLOW Closed TUR irrigation bag system UB 4000TUR with pump ball for controlled post-operative bladder irrigation

Capacity: 4000 ml

Product specifications:

- ✓ Can be fixed safely and correctly to many different types of bed frames by means of the twin holder with fixing cord.
- ✓ Bacteria- and fluid-tight venting filters ensure permanent pressure compensation and provide for optimal drainage in the system.
- ✓ Due to controlled irrigation of the bladder with the pump ball integrated in the feeding tube, post-operative tamponade of the bladder can be avoided.
- ✓ A doghouse on the bag ensures safe and hygienic positioning of the drain tap.



Capacity	Tube length	Article number	PZN
4000 ml	120 cm	04/02/4000/TUR/GHC	06189004



Pump ball

The specific interaction of the pump ball and the two tube clamps allows to flush out blood clots in the bladder and tube system.



Profiled tube

The octagonal inner profile of the feeding tube assures excellent stability, but at the same time it is also flexible. A high lumen guarantees a safe drainage into the drip chamber in any situation.



Drip chamber

An anti-reflux valve divides the drip chamber and reliably prevents a backflow of urine and thereby obviates the risk against an ascending infection.



Drain tap in the doghouse

The special design of the piston in the drain tap ensures safe and easy operation and prevents leakages caused by improper use of the piston.

Packaging Unit: 10 pieces

Packaging: sterile

Pictures may differ from the original