

# CARE FLOW UB 1000 Closed urine drainage systems

## Leg bag systems for day & night

Capacity: 1000 ml

### Product specifications:

- ✓ Closed urine drainage system
- ✓ Bacteria- and fluid-tight venting filters provide for permanent pressure compensation and optimal drainage in the system.
- ✓ A doghouse on the bag ensures safe and hygienic positioning of the drain tap.
- ✓ T-drain tap.
- ✓ Needle free sampling port.
- ✓ Clear gradvation marks.



Capacity	Tube length	Article number	PZN
1000 ml	30 cm	04/03/1000/30/GHC	06718885
1000 ml	60 cm	04/03/1000/60/GHC	06718891
1000 ml	90 cm	04/03/1000/90/GHC	06718916



#### Stepped connector

The special design of the connector steps provides for a reliable connection with the ring profile in the funnel end of the catheter. Needle-free urine sampling avoids needle injuries.



#### Fleece coating

The backside of the leg bag is coated with sweat-absorbing fleece. This effectively prevents irritation of the skin and considerably improves user comfort.



#### Drip chamber

An anti-reflux valve in the drip chamber provides for reliable prevention of urine reflux and highly efficient protection against ascending migration of bacteria.



#### Drain tap

The special design of the piston in the drain tap ensures safe and easy operation and prevents leak-ages caused by improper use of the piston.

**Packaging Unit:** 30 pieces

**Packaging:** sterile

**Fixation for leg bags on page 51-56.**

Pictures may differ from the original