

CARE FLOW UB 2000E, UB 2000US

Sterile / non-sterile single-use bedside bag with holder

Capacity: 2000 ml

Product specifications:

- ✓ The system is equipped with a universal staggered cone and can be connected to all common catheters.
- ✓ The product must be changed daily.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.
- ✓ Use the twin holder shown below to attach the bag to many different types of bed frames.
- ✓ The urine bag is emptied safely by means of the T-valve.



Article	Capacity	Sterility	Tube length	Article number	PZN
UB 2000E	2000 ml	sterile	90 cm	04/04/2000/90/GHC	16047274
UB 2000US	2000 ml	non-sterile	90 cm	04/06/2000/90/GHC	12418549
UB 2000US	2000 ml	non-sterile	130 cm	04/06/2000/130/GHC	08845487

45

UB 2000E: **Packaging Unit:** 100 pieces, **Packaging:** sterile
 UB 2000US: **Packaging Unit:** 100 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces

Holder for UB 2000E and UB 2000US

- ✓ Plastic hanger for bedside bags with twin holder. Equipped with two pins taking up single-use bedside bags hanger holes.



Article number	PZN
04/99/BHL/PVC/GHC	06718951

Packaging Unit: 100 pieces **Sub-packaging:** 10 pieces each
Packaging: non-sterile

Pictures may differ from the original