

# CARE FLOW UB 2000A Disposable bedside bag

## Sterile bedside bag with integrated twin holder

Capacity: 2000 ml

### Product specifications:

- ✓ Can be fixed safely and correctly to many different types of bed frames by means of the twin holder.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.
- ✓ Clear graduation markings allow for reliable volume control.
- ✓ A slit ensures the safe and hygienic positioning of the dRAINTAP.
- ✓ T-drain tap.
- ✓ Kink resistant tube.



| Capacity | Tube length | Article number     | PZN      |
|----------|-------------|--------------------|----------|
| 2000 ml  | 110 cm      | 04/04/2000/110/GHC | 08871289 |

44



#### Staggered connector

The staggered connector provides a reliable connection with the catheter. A needle-free sampling of urine is possible at any time, due to the needle-free sample port, which prevents needlestick injuries.



#### Profiled tube

The octagonal inner profile of the feeding tube assures excellent stability, but at the same time it is also flexible. A high lumen guarantees a safe drainage into the drip chamber in any situation.



#### Venting filter

A bacteria and liquid-tight venting filters ensure permanent pressure equalisation. This ensures an optimal drainage of the system.



#### Drain tap

The special design of the piston in the drain tap ensures safe and easy operation and prevents leakages caused by improper use of the piston.

**Packaging Unit:** 50 pieces  
**Packaging:** sterile

Pictures may differ from the original