

CARE FLOW UB 2000USo, UB 1500USo, SB 1500USo

Non-sterile single-use bedside bag/secretion bag without drain tap

Capacity: 2000 ml, 1500 ml

Product specifications:

- ✓ The system is equipped with a universal staggered cone and can be connected to all common catheters.
- ✓ The product must be changed daily.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria. (except the SB 1500USo)
- ✓ Use the twin holder shown below to attach the bag to many different types of bed frames.



Article	Capacity	Tube length	Anti-reflux valve	Article number	PZN
UB 2000USo	2000 ml	90 cm	yes	04/07/2000/90/GHC	08845501
UB 1500USo	1500 ml	90 cm	yes	04/07/1500/90/01	16839791
SB 1500USo	1500 ml	90 cm	no	04/07/1500/90/00	16839785

46

UB 2000USo: **Packaging Unit:** 150 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces
 UB 1500USo: **Packaging Unit:** 200 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces
 SB 1500USo: **Packaging Unit:** 200 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces

Holder for UB 2000USo, UB 1500USo and SB 1500USo

- ✓ Plastic hanger for bedside bags with twin holder. Equipped with two pins taking up single-use bedside bags hanger holes.

Article number	PZN
04/99/BHL/PVC/GHC	06718951

Packaging Unit: 100 pieces **Sub-packaging:** pieces each
Packaging: non-sterile



Pictures may differ from the original