

CARE FIX

Velcro strap premium set

Length can be individually adjusted

Product specifications:

- ✓ Sets consists of:
an extra wide holding band with two sewn-on Velcro straps to fixate the upper part of the leg bag.
- ✓ Both sets are equipped sewn-on loop to fixate the leg bag tube.
- ✓ The second Velcro Strap ensures a safe fixation of the lower part of the leg bag.
- ✓ Silicone strips on the inner side of the holding band ensures a good grip on the leg.
- ✓ The elastic material prevents a disruption of the blood circulation on the thigh or lower leg.
- ✓ A 70 cm long and individually adjustable Velcro strap can be worn on both the thigh as well as the lower leg.
- ✓ The attached loop for an easier opening of the holding band.



Lenght	Article number	PZN
70 cm	05/01/BKL/PRE/SET	09241099

Packaging Unit: 10 pieces

Packaging: non-sterile

52

- ✓ This Velcro strap premium set can be used for the following urine leg bag systems: UB1000, UB750/UB750US, UB750F/UB750UF, UB500/UB500US.
- ✓ The premium set with buttons (Art.-Nr. 05/01/BKL/PRE/03) comes with additional four buttons to fixate the upper part of the leg bag.

Care instructions

- 60°
- Do not put in the dryer
- Do not iron

Lenght	Article number	PZN
70 cm	05/01/BKL/PRE/03	12357931

Packaging Unit: 10 pieces

Packaging: non-sterile



Pictures may differ from the original