

CARE FIX

Universal Velcro strap




Length can be individually adjusted

Product specifications:

- ✓ Universal leg bag holding band with two attached Velcro straps to go through the upper eyelets of the leg bag and a special sewn on fixation for the tube of the leg bag system.
- ✓ The elastic material ensures that disturbances of the blood circulation in the thigh and lower leg are prevented as far as possible.
- ✓ The Velcro strap is 70 cm long and can be cut to the individual length. So it can be used on the thigh and lower leg.
- ✓ This universal Velcro strap can be used for the following urine leg bag systems: UB1000, UB750/UB750US, UB750F/UB750UF, UB500/UB500US.



Care instructions:

-  60°
-  Do not put in the dryer
-  Do not iron



Length	Article number	PZN
70 cm	05/01/BKL/Uni/GHC	06718968

Packaging Unit: 30 pieces
Packaging: sterile

Pictures may differ from the original