

CARE FIX

Thigh cuffs sizes: XL, L, M

Product specifications:

- ✓ The fabric leg bag holder made from a washable cotton blend has an integrated leg bag pocket.
- ✓ A silicon strip on the inner side ensures a good grip.
- ✓ The fabric leg bag holder can be worn on both the right as well as the left thigh.
- ✓ The bags are positioned into the integrated pocket with a Velcro fastener; including an opening for the drain tap.
- ✓ The CARE FIX cuffs can be used for the following urine leg bag systems:
UB1000, UB750/UB750US, UB750F/
UB750UF, UB500/UB500US.



54

Care instructions:

- 60°
- Do not put in the dryer
- Do not iron

Oberschenkelstulpen

Size	Article number	PZN
XL	05/02/12/XL/GHC	09080408
L	05/02/12/L/GHC	09080383
M	05/02/12/M/GHC	09080377

Measurement: XL = 61 cm, L = 58 cm, M = 54 cm

Packaging Unit: 10 pieces

Packaging: non-sterile



Pictures may differ from the original