

CARE FIX

Lower leg cuffs sizes: L, M, S

Product specifications:

- ✓ The fabric lower leg bag holder made from a washable cotton blend has an integrated leg bag pocket.
- ✓ A silicon strip on the inner side ensures a good grip.
- ✓ The fabric lower leg bag holder can be worn on both the right as well as the left lower leg.
- ✓ Leg bags are positioned into the integrated pocket with a Velcro fastener; including an opening for the drain tap.
- ✓ The CARE FIX cuffs can be used for the following urine leg bag systems: UB1000, UB750/UB750US, UB750F/UB750UF, UB500/UB500US.



Care instructions:

- 60°
- Do not put in the dryer
- Do not iron

Unterschenkelstulpen

Size	Article number	PZN
L	05/02/13/L/GHC	09080437
M	05/02/13/M/GHC	09080420
S	05/02/13/S/GHC	09080414

Measurement: L = 37 cm, M = 35 cm, S = 32 cm

Packaging Unit: 10 pieces

Packaging: non-sterile

