

Curan Complete & Closed System

IC disposable catheter with integrated collection bag, for women & men




IC disposable catheter

- ✓ Complete systems with disposable catheter and integrated collection bag.
- ✓ Good gliding properties thanks to hydrogel coating.
- ✓ Polished drainage eyes enable comfortable insertion and prevent irritation or injury to the urethral mucosa.
- ✓ Disposable catheters are ready to use.
- ✓ The catheter material is phthalate-free.
- ✓ The introducer tip can be inserted hygienically into the urethra to bypass the germ-colonised urethral orifice by not coming into contact with the catheter itself.



Curan Complete compact system
Collection bag: 1.000 ml



Curan Closed System
Collection bag: 1.500 ml

Curan Complete with collection bag, 1.000ml

Size	Lenght	Article number	PZN
CH 12	40 cm	CURAN-CC12	19341933
CH 14	40 cm	CURAN-CC14	19341956

Packaging Unit: 30 pieces

Packaging: sterile

Curan Closed System with collection bag, 1.500ml

Size	Lenght	Article number	PZN
CH 10	40 cm	CURAN-CB10	19341850
CH 12	40 cm	CURAN-CB12	19341867
CH 14	40 cm	CURAN-CB14	19341896
CH 16	40 cm	CURAN-CB16	19341904

Packaging Unit: 25 pieces

Packaging: sterile



Insertion tip