

CARE CATH_{Leiv}

IC disposable catheter, hydrophilic, 41 cm
Ready for immediate use, Nelaton

Product specifications:

- ✓ Very good gliding properties thanks to hydrophilic coating.
- ✓ Ready for immediate use thanks to ready-to-use, hydrophilic coating.
- ✓ Atraumatic application due to rounded drainage eyes.
- ✓ Non-contact application and good guideability due to the correct rigidity of the catheter.
- ✓ Good grip of the handle allows good handling, even with limited dexterity.
- ✓ Easy connection of urine collection bag and/or instillation medication via the handle cone.



Atraumatic, rounded drainage eyes




CARE CATH_{Leiv}

Size	Length	Item	number
CH 12	41 cm	LEIV-H1241	18890833
CH 14	41 cm	LEIV-H1441	18890856
CH 16	41 cm	LEIV-H1641	18890862

Packaging unit: 30 pieces

Packaging: sterile



Grip of the handle

Illustrations may differ from the original