

Product specifications:

- ✓ Very good gliding properties thanks to hydrophilic coating.
- ✓ Ready for immediate use thanks to ready-to-use, hydrophilic coating.
- ✓ Atraumatic application due to rounded drainage eyes.
- ✓ Optimal length of 17 cm ensures complete bladder emptying.
- ✓ Non-contact application and good guideability due to the correct rigidity of the catheter.
- ✓ Good grip of the handle allows good handling, even with limited dexterity.
- ✓ Easy connection of urine collection bag and/or instillation medication via the handle cone.



Atraumatic, rounded drainage eyes



Size	Length	Item	number
CH 10	17 cm	LIV-H1017	18890885
CH 12	17 cm	LIV-H1217	18890916
CH 14	17 cm	LIV-H1417	18890922

Packaging unit: 30 pieces

Packaging: sterile



Grip of the handle

Illustrations may differ from the original